

**MIBI 20X TBS-T**

Catalog: 567005

Storage: Store at 4°C

**Description** Tris buffered saline (TBS), pH 7.6, with the detergent Tween® 20 is intended to be used as a wash buffer to remove nonspecific antibody binding and as a blocking buffer diluent in the MIBI Staining Protocol.

**Directions For Use** Dilute 20X TBS to a 1X solution using MIBI diH<sub>2</sub>O.

**Validation** Each lot is quality control tested to ensure the absence of contaminating metal impurities.